# Vacuum furnace operation manual

Juhee Ko

Department of Mechanical Engineering  
Korea Advanced Institute of Science and Technology

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

# Panel configuration

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ACE-VAC VACUUM FURNACE

The top section of the panel features three digital control units. On the left is a 'PROGRAM CONTROLLER' (YOKOGAWA PV) with a green display showing '55P'. In the center is an 'OVER TEMPERATURE' (KX3N) unit with red displays showing '88.14' and '1250'. On the right is a 'VACUUM GAUGE' (KVC5000N) with green displays showing '76.2' and 'Err'. Below these units are the labels '과열 방지' (Overheat prevention) and '압력 게이지' (Pressure gauge).

The middle section contains a row of control elements. From left to right: 'ROUGH PUMP ON/OFF' (green indicator), 'DIFFUSION PUMP ON/OFF' (black knob), 'HEATER OFF/ON' (black knob), 'ROUGH VALVE CLOSE/OPEN' (black knob), 'FORE VALVE CLOSE/OPEN' (green indicator), and 'HI-VAC VALVE CLOSE/OPEN' (black knob). Below these are six red indicator lights labeled 'WATER LOW', 'AIR LOW', 'OVER TEMP', 'SCR FAULT', 'BUZZER', and 'B/Z STOP'. The label '펌프 스위치' (Pump switch) is centered above the top row, and '각종 경고등' (Various warning lights) is centered below the bottom row.


The bottom section of the panel features four analog meters. From left to right: 'MAIN V-METER' (0-200V), 'MAIN A-METER' (0-60A), 'HEATER V-METER' (0-50V), and 'HEATER A-METER' (0-1500A).

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

Wait 5 s

Wait until pressure drops under  $5 * 10^{-1}$  Torr (~ 10 min) and turn off the valve

Wait 5 s

Wait until target pressure is reached

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(RP 30 min, DP 30 min)

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# Temperature controller operation (programming)

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